




Reusable Laparoscopic Instruments



Trocar

Product	Description	Dimension	Code
	Pyramidal Trocar	φ5X95	221-51155
		φ10X95	221-01155
		φ12X95	221-81155
		φ5X150	221-53155
		φ10X150	221-03155
		φ12X150	221-83155
	Auto Shield Trocar	φ5X95	221-51255
		φ10X95	221-01255
		φ12X95	221-81255
		φ5X150	221-53255
		φ10X150	221-03255
		φ12X150	221-83255
	Pyramidal Trocar with Screw	φ5X95	221-51165
		φ10X95	221-01165
		φ12X95	221-81165
	Auto Shield Trocar with Screw	φ5X95	221-51265
		φ10X95	221-01265
		φ12X95	221-81265
	Hasson Blunt Tip Trocar	φ10X95	221-01555

Trocar



Product	Description	Dimension	Code
	Veress Needle	φ2X120	221-21127
		φ2.7X150	221-27157
	Reducer Sleeve	φ10-φ5	221-01785
		φ12-φ5	221-81785
	Trocar Stop	φ5X35	221-50915
		φ10X35	221-00915
		φ12X35	221-80915
	Sealing Cap	φ5	221-50925
		φ10	221-00925
		φ12	221-80925
	Reducer	φ10-φ5	221-00725
		φ12-φ5	221-80725
	Reducer	φ10-φ5	221-00715

Reusable Laparoscopic Instruments



Safe

Insulated sleeve exposes a shorter metal base to avoid the accidental damage caused by electrical leakage.



Durable

Strictly selected materials plus special production processes greatly extend the product's life.



Ergonomic

One-click to disassemble into three parts for more thorough cleaning and sterilization.



Dissecting and Grasping Forceps

- **Diamond serration clamp teeth**
Stronger grasping force with less tissue damage
- **No obvious bulge at the opening joint**
Reduce blood and tissue residues, easy to clean and reduce cross infection
- **Shorter metal base**
Less unintentional electrical injury

Product	Description	Dimension	Code
	KELLY Dissecting and Grasping forceps, Maryland	φ5X360	222-5S3013
		φ5X360, With Ratchet	222-5T3013
		φ5X430	222-5S4013
		φ5X430, With Ratchet	222-5T4013
	KELLY Dissecting and Grasping forceps, long, Maryland	φ5X360	222-5S3023
		φ5X360, With Ratchet	222-5T3023
		φ5X430, With Ratchet	222-5T4023
	Dissecting and Grasping forceps, right-angled	φ5X360	222-5S3123
		φ5X360, With Ratchet	222-5T3123
		φ5X430, With Ratchet	222-5T4123
	Straight Dissecting and Grasping forceps, atraumatic	φ5X360	222-5S3043
		φ5X360, With Ratchet	222-5T3043
	Dissecting and Grasping forceps, atraumatic	φ10X360	222-0S3103
		φ10X360, With Ratchet	222-0T3103
	Dissecting and Grasping forceps, right-angled	φ10X360	222-0S3133
		φ10X360, With Ratchet	222-0T3133

Grasping Forceps



No Ratchet



With Ratchet

Product	Description	Dimension	Code
	Bowel Grasper, fenestrated, short	φ5X360	222-5S3213
		φ5X360, With Ratchet	222-5T3213
	Atraumatic Grasping forceps	φ5X360	222-5S3633
		φ5X360, With Ratchet	222-5T3633
		φ5X430, With Ratchet	222-5T4633
	CROCE-OLMI Grasping forceps, atraumatic, curved	φ5X360	222-5S3453
		φ5X360, With Ratchet	222-5T3453
		φ5X430, With Ratchet	222-5T4453
	Grasping forceps, w specially fine atraumatic serration	φ5X360	222-5S3473
		φ5X360, With Ratchet	222-5T3473
		φ5X430	222-5S4473
		φ5X430, With Ratchet	222-5T4473
	Bowel Grasper, fenestrated	φ5X360	222-5S3303
		φ5X360, With Ratchet	222-5T3303
	BABCOCK Grasping forceps	φ5X360	222-5S3553
		φ5X360, With Ratchet	222-5T3553
	Grasping forceps, 2X4 teeth	φ5X360	222-5S3393
		φ5X360, With Ratchet	222-5T3393

Grasping Forceps



No Ratchet



With Ratchet

Product	Description	Dimension	Code
	Grasping forceps, atraumatic	φ5X360	222-5S3333
		φ5X360, With Ratchet	222-5T3333
	VANCAILLIE Adhesion Forceps, one jaw fenestrated	φ5X360	222-5S3493
		φ5X360, With Ratchet	222-5T3493
	MANHES biopsy forceps	φ5X360	222-5S3903
		φ5X360, With Ratchet	222-5T3903
	MANHES Grasping forceps, "tiger-jaws", 2X4 teeth	φ5X360	222-5S3363
		φ5X360, With Ratchet	222-5T3363
	Claw forceps, 2X3 teeth, short	φ10X360	222-0S3403
		φ10X360, With Ratchet	222-0T3403
	SAWALHE Tissue Grasping forceps	φ10X360	222-0S3423
		φ10X360, With Ratchet	222-0T3423
	BABCOCK Grasping forceps	φ10X360	222-0S3573
		φ10X360, With Ratchet	222-0T3573

Grasping Forceps



No Ratchet



With Ratchet

Product	Description	Dimension	Code
	Claw forceps, 2X3 teeth	φ10X360	222-0S3413
		φ10X360, With Ratchet	222-0T3413
	Spoon forceps	φ10X360	222-0S4883
		φ10X360, With Ratchet	222-0T4883
	VANCAILLIE Oviduct Forceps	φ5X360	222-5S3583
		φ5X360, With Ratchet	222-5T3583
	Dissecting and Grasping forceps, "alligator jaws"	φ5X360	222-5S3773
		φ5X360, With Ratchet	222-5T3773
	DeBAKEY Grasping forceps, curved and slender jaws	φ5X360	222-5S3813
		φ5X360, With Ratchet	222-5T3813
	Grasping forceps, atraumatic, spoon-shaped	φ5X360	222-5S3223
		φ5X360, With Ratchet	222-5T3223
	MANHES Dissecting and Grasping forceps, "duckbill"	φ5X360	222-5S3283
		φ5X360, With Ratchet	222-5T3283


Monopolar Scissors




- Sharp**
 The sharpest blade enables fast cutting.
- Smooth**
 Involute design of blade line ensures more smooth cutting.
- Durable**
 The vacuum heat treatment process makes sure that the metal hardness control accuracy reaches HRC1 (Rockwell hardness), increasing the metal life by 3-4 times compared with the traditional method.

Product	Description	Dimension	Code
	METZENBAUM Scissors, curved	φ5X360	222-5S3022
		φ5X430	222-5S4022
	Scissors, spoon-shaped blades, serrated, curved	φ5X360	222-5S3012
		φ5X430	222-5S4012
	Scissors, straight	φ5X360	222-5S3092
		φ5X430	222-5S4092

Bipolar Dissecting Forceps



Integrating dissecting, grasping, and coagulating together helps improve surgery efficiency by reducing the frequency of instruments changes.







- Strong grasping:**
 Enhanced force with wider opening angle to achieve better grasping.
- Guaranteed Safety:**
 Advanced insulation design guarantees the safety and reliability of bipolar coagulation.
- Innovative anti-sticking:**
 Sticking prevention with coated jaw.

Product	Description	Dimension	Code
	Bipolar Straight Dissecting Forceps	φ5X360	222-5S3011
		φ5X430	222-5S4011
	Bipolar Curved Maryland Dissecting Forceps	φ5X360	222-5S3021
		φ5X430	222-5S4021



Monopolar Electrodes





Product	Description	Dimension	Code
	Hook Tip Coagulation Electrode	φ5X330	222-51110
	Ball Tip Coagulation Electrode	φ5X330	222-51220
	Spatula Tip Coagulation Electrode	φ5X330	222-51440
	Rod Tip Coagulation Electrode	φ5X330	222-51610

Cables



All cables are autoclavable.

Product	Description	Dimension	Code
	Monopolar cable	3200mm	999-030
	Bipolar Cable	3200mm	999-031


ENDO-Retractor

Product	Description	Dimension	Code
	ENDO-Retractor	φ5X360	221-53013
	ENDO-Retractor	φ5X360	221-53023


Fans Dissecting Forceps

Product	Description	Dimension	Code
	3 Finger Retractor	φ5X330	221-55923
	5 Finger Retractor	φ10X330	221-05923

Myoma Drill

Product	Description	Dimension	Code
	Myoma Drill	φ5X330	221-53109
		φ5X420	221-54109
		φ10X330	221-03109
		φ10X420	221-04109

Clip Appliers




Reliable ligation of blood vessels to avoid intraoperative and postoperative bleeding.

Product	Description	Dimension	Code
	Endo Applier for Medium Large Polymer Clip	φ5X330	221-55636
	Endo Applier for Large Polymer Clip	φ10X330	221-05646
	Endo Applier for Extra Large Polymer Clip	φ10X330	221-05656





Clip Appliers

Product	Description	Dimension	Code
	Endo Applier for Small Titanium Clip, Single action	φ5X330	221-55116
	Endo Applier for Medium Titanium Clip, Single action	φ10X330	221-05126
	Endo Applier for Medium Large Titanium Clip, Single action	φ10X330	221-05136
	Endo Applier for Large Titanium Clip, Single action	φ10X330	221-05146
	Open Surgery Applier for Small Titanium Clip	28 cm	221-00516
	Open Surgery Applier for Small Titanium Clip	18 cm	221-01516
	Open Surgery Applier for Medium Titanium Clip	28 cm	221-00526
	Open Surgery Applier for Medium Titanium Clip	23 cm	221-01526
	Open Surgery Applier for Medium Large Titanium Clip	28 cm	221-00536
	Open Surgery Applier for Large Titanium Clip	28 cm	221-00546

Needle Holders




Needle Holders

Product	Description	Dimension	Code
	V Handle Straight Needle Holder	φ5X330	221-56416
	V Handle Left Curved Needle Holder	φ5X330	221-56426
	Pistol Handle Straight Needle Holder	φ5X330	221-56516
	Pistol Handle Left Curved Needle Holder	φ5X330	221-56526
	Pistol Handle Right Curved Needle Holder	φ5X330	221-56536
	Pistol Handle Needle Holder with Reposition Function	φ5X330	221-56546
	O Handle Straight Needle Holder	φ5X330	221-56816
	O Handle Left Curved Needle Holder	φ5X450	221-56876
	O Handle Right Curved Needle Holder	φ5X330	221-56826
	O Handle Needle Holder with Reposition Function	φ5X450	221-56886
	O Handle Needle Holder with Reposition Function	φ5X330	221-56836

Suture Fascia Fix Forceps

Product	Description	Dimension	Code
	Suture Fascia Fix Forceps	φ3X200	221-32013


Knot Pusher

Product	Description	Dimension	Code
	Knot Pusher	φ4X330	221-43019
		φ4X330	221-43029

Suction & Irrigation Apparatus






Three ergonomic handles satisfy different customers' preference.



Strong suction and irrigation, as well as blunt dissection fulfill different operation demands.

Suction & Irrigation Apparatus

Product	Description	Dimension	Code
	Suction and irrigation	φ5X330	221-53114
	Suction and irrigation	φ5X330	221-53224
		φ10X330	221-03224
	Suction and irrigation	φ5X330	221-53334
		φ10X330	221-03334

Puncture Injection Needle



Description	Dimension	Code
Puncture Injection Needle (7#)	φ5X330	221-53015
Puncture Injection Needle (9#)	φ5X330	221-53025
Puncture Injection Needle (12#)	φ5X330	221-53035
Puncture Injection Needle (14#)	φ5X330	221-53045
Puncture Injection Needle (15#)	φ5X330	221-53055
Puncture Injection Needle (18#)	φ5X330	221-53065